

E. coli O157 Antiserum



pro-lab.com/resources



60

REF PL.6250

INSTRUCTIONS FOR USE

INTENDED USE

Pro-Lab *E. coli* O157 Antiserum is for use in the slide agglutination test for the presumptive identification of *Escherichia coli* serotype O157 antigen on laboratory culture media.

SUMMARY AND EXPLANATION

Escherichia coli serotype O157:H7 is a verotoxin producing (VT- producing) pathogen.^{1,2} This serotype has been reported as an etiological agent in sporadic and outbreak cases of haemorrhagic colitis.^{3,4,5} It is also associated with haemolytic uraemic syndrome.⁶ Certain *E. coli* serotypes other than O157:H7 also produce verotoxin.^{7,8,9} However, the diarrhoea caused by these other serotypes is not usually bloody. Additionally, *E. coli* serotype O157:H7 does not ferment sorbitol whereas the majority of other serotypes do ferment sorbitol.^{10,11} Therefore, if Sorbitol-MacConkey agar medium is used as a primary screen, the colonies of *E. coli* serotype O157:H7 appear colourless (non-sorbitol fermenting colonies - NSFC) while colonies of other serotypes appear characteristically pink (sorbitol fermenting colonies - SFC).¹¹

The work of Kauffmann¹², Edward and Ewing¹³, Ewing¹⁴ and Orskov¹⁵ contributed to the development of a system for serological typing of *E. coli* cultures and resulted in an antigenic classification scheme which can be used to identify the serotypes of *Escherichia coli* which are associated with bacteriuria or diarrheal disease.

The principle of the test involves mixing the suspected organisms with the antiserum containing *E. coli* O157 antibodies. The bacteria will agglutinate (clump) in the presence of homologous antiserum.

PRINCIPLE OF THE TEST

The principle of the test involves mixing the suspected organisms with the antiserum containing *E. coli* O157 antibodies. The bacteria will agglutinate (clump) in the presence of homologous antiserum.

MATERIALS PROVIDED

Pro-Lab *E. coli* O157 Antiserum is prepared using delipidized, whole absorbed rabbit serum containing antibodies to *E. coli* serotype O157.

Pro-Lab *E. coli* O157 Antiserum is supplied in a dropper bottle containing 3.0 ml of ready-to-use diluted antiserum with 0.01% thimerosal as preservative.

MATERIALS REQUIRED BUT NOT PROVIDED

- Glass Slides
- Normal saline (0.85% sodium chloride) or Phosphate Buffered Saline (PBS)
- Inoculating loop or needle
- Mixing sticks

STABILITY AND STORAGE

Pro-Lab *E. coli* O157 Antiserum should be stored tightly capped at 2-8°C. Stored under these conditions will be stable until the expiry date shown on the product label.

WARNINGS AND PRECAUTIONS

1. For in vitro diagnostics use.
2. Do not use the antiserum after the expiry date shown on the product label.
3. The antiserum contains thimerosal, which is a highly toxic mercury based compound. Although the amount of thimerosal in the antiserum is minimal, safety precautions should be taken in handling, processing and discarding the reagent.
4. Avoid contamination of the reagent bottle.
5. Universal precautions should be taken in handling, processing and discarding all clinical specimens. All test materials should be considered potentially infectious during and after use and should be handled and disposed of appropriately.
6. The procedures, storage conditions, precautions and limitations specified in these directions must be adhered to in order to obtain valid test results.
7. Antiserum contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

SPECIMEN COLLECTION AND PREPARATION

Clinical specimens should be cultured on Sorbitol-MacConkey medium. NSFC may be subcultured on non-selective agar medium. Colonies from overnight growth must be cleanly removed from agar surface for testing using a sterile loop. Young, fast growing cultures will yield typical results.

TEST PROCEDURE

1. Place two separate drops of normal saline or phosphate buffered saline on a clean glass slide.
2. Take a suspect *Escherichia coli* colony from an overnight culture plate and mix thoroughly with both drops of normal saline or phosphate buffered saline on the slide to obtain a smooth suspension.
3. Add one drop of antisera to one of the bacterial suspension drops on the slide, to the other (control) add one drop of normal saline or phosphate buffered saline.
4. Mix the antiserum with the bacterial suspension using a mixing stick. Then mix the saline or phosphate buffered saline (control) with a fresh mixing stick.
5. Gently rock the slide back and forth for one minute and observe for agglutination under normal lighting conditions or using a low power objective.

INTERPRETATION OF RESULTS

A distinct agglutination (granular clumping) in the antiserum test, within 60 seconds, is regarded as a positive result. There must be no agglutination in the saline or phosphate buffered saline control or else the test is not valid (autoagglutination).

LIMITATIONS OF THE PROCEDURE

1. A normal saline or phosphate buffered saline control should be included in every test to insure the specificity of the reaction.
2. Rough strains give autoagglutination in slide tests. False positives usually agglutinate in control saline or phosphate buffered saline.
3. It is recommended to check the potency of *Escherichia coli* antisera with stock cultures of known antigenic structure.
4. The antiserum is a presumptive identification or confirmation of cultures which have been previously characterized biochemically.

AVAILABILITY












Cat. No	Description
PL.070B	Prolex <i>E. coli</i> O157 Latex Test Reagent Kit (50 Tests)
PL.071B	Prolex <i>E. coli</i> O157 Latex Test Reagent Kit (100 Tests)

REFERENCES

1. Konowalchuk J., Speirs J.I., Stavric S. 1977. Vero response to a cytotoxin of *Escherichia coli*. Infect. Immun. 18:775-779.
2. Ratnam S., March S.B., Ahmed R., Bezanson G.S., Kasatiya S. 1988. Characterization of *Escherichia coli* serotype O157:H7. J. Clin. Microbiol. 26:2006-2012.
3. C.D.C. 1982. Isolation of *E. coli* O157:H7 from sporadic cases of hemorrhagic colitis. United States MMRW 31:580-585.

4. Johnson W.M., Lior H., Bezanson G.S. 1983. Cytotoxic *Escherichia coli* 0157:H7 associated with haemorrhagic colitis in Canada. Lancet i:76.
5. Krishnan C., Fitzgerald V., Dakin S., Behme R.J. 1987. Laboratory investigation of outbreak of hemorrhagic colitis caused by *Escherichia coli* 0157:H7. J. Clin. Microbiol. 25:1043-1047.
6. Karmali M.A., Steele B.T., Petric M., Lim C. 1983. Sporadic cases of haemolytic-uraemic syndrome associated with faecal cytotoxin and cytotoxin-producing *Escherichia coli* in stools. Lancet. i:619- 620.
7. Karmali M.A., Petric M., Lim C., Cheung R., Arbus G.S. 1985. Sensitive method for detecting low numbers of verotoxin-producing *Escherichia coli* in mixed cultures by use of colony sweeps and polymyxin extraction of verotoxin. J. Clin. Microbiol. 22:614-619.
8. Law D. 1988. Virulence factors of enteropathogenic *Escherichia coli*. J. Med. Microbiol. 26:1-10.
9. Scotland S.M., Day N.P., Rowe B. 1980. Production of a cytotoxin affecting vero cells by strains *Escherichia coli* belonging to traditional enteropathogenic serogroups. FEMS Microbiol. Lett. 7:15- 17.
10. Farmer III J.J., Davis B.R. 1985. H7 Antiserum-sorbitol fermentation medium: a single tube screening medium for detecting *Escherichia coli* 0157:H7 associated with hemorrhagic colitis. J. Clin. Microbiol. 22:620-625.
11. March S.B., Ratnam S. 1986. Sorbitol-MacConkey medium for detection of *Escherichia coli* 0157:H7 associated with hemorrhagic colitis. J. Clin Microbiol. 23:869-872.
12. Kauffmann, F. 1947. J. Immunology 57:71-100.
13. Edwards, P.R. and Ewing, W.H. 1972. Identification of Enterobacteriaceae. 3rd edition. Burgess. Minneapolis, Minnesota.
14. Ewing, W.H. 1969. Public Health Lab. 27:19-30.
15. Orskov, F. 1956. Acta. Pathol. Microbiol. Scand. 29:373.

SYMBOLS GLOSSARY

Symbol	Meaning
	Manufacturer
	Use-by date
	Lot number
	Catalogue number
	Temperature limit
	Consult instructions for use or consult electronic instructions for use
	In vitro diagnostic medical device
	Contains sufficient for <n> tests
	Indicates European Conformity
	Indicates UK Conformity
	Authorized representative in the European Community / European Union

Pro-Lab Diagnostics U.K.

3 Bassendale Road, Bromborough,
Wirral, Merseyside CH62 3QL
United Kingdom
Tel.: 0151 353 1613
www.pro-lab.co.uk
E-mail: uksupport@pro-lab.co.uk




Pro-Lab Diagnostics Canada

20 Mural Street, Unit #4
Richmond Hill, Ontario L4B 1K3
Canada
Tel.: 905-731-0300
www.pro-lab.com
E-mail: support@pro-lab.com

Pro-Lab Diagnostics U.S.A.

1301 Blue Ridge Drive, Suite #101
Georgetown, Texas 78626-3269
United States
Tel.: 512-832-9145
www.pro-lab-direct.com
E-mail: support@pro-lab.us

 Advena Ltd. Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013 Malta